



## Which Z1 Is Right For You?

Z1 CPAP

Z1 Auto

Fixed Pressure Mode	●	●
Auto Adjusting Pressure Mode		●
Auto Stop Feature	●	
Auto Start/Stop Feature		●
Approved for In-Flight Use	●	●
PC Data Capable	●	●
Mobile Data Capable		●
Pressure Range 4-20 cm H <sub>2</sub> O	●	●
Ramp 0-45 minutes	●	●
Z-Breathe® Comfort Feature	●	●
Qtube Included*	●	●
No Special Mask Required	●	●
Auto Adjusting Voltage for World Travel (100-240v)	●	●
Weight: Z1 Auto: 299g, Z1 CPAP: 292g		
Dimensions: 16.5x8.4x5.1cm		

© 2019 Breas Medical – All rights reserved. Breas Medical reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your Breas representative for the most current information. Breas and the Breas logo are trademarks of Breas Medical AB. COR-SAL-MAR-001112-Rev5

# BREAS

Breas Medical AB · Företagsvägen 1 · SE-435 33 Mölnlycke · Sweden  
Phone +46 31 86 88 00 · [breas.com](http://breas.com)

**The PowerShell® - Totally Integrated Battery Power to-go.** The optional PowerShell® accessory is a completely integrated cpap battery solution that delivers a full night's sleep on a single charge. No cords required.



### Battery Specifications

Battery capacity	99.4 Wh
Operating temperature	5°-35°C (41°-95°F)

### Modular Batteries


Need more nights of off-the-grid power? Buy an additional battery module. Boom. Have fun.



## Nitelog™ Mobile App

Exclusively for the Z1™ Auto

Interested in your sleep data? Download the Nitelog™ app from the iTunes App Store® or Google™ play and have all of your sleep data available at your fingertips. Nitelog uses the Z1™ Auto's built-in Bluetooth chip to sync your data wirelessly so it is very simple to set up and use. Nitelog also enables sharing by email so that you can easily send your data to your doctor.



# Z1™ auto

Ultra-Small, Travel CPAP.





Z1<sup>TM</sup> auto

# Take Your Adventures Off the Grid

Whether you need the fixed pressure Z1 CPAP or the auto-adjusting Z1 Auto, you can take your CPAP anywhere you want to go.

The optional PowerShell® with Extended Life Battery provides battery power for a full night's sleep on a single charge. It's great for camping or any where access to power is limited.

## Travel Anywhere in the World

Both Z1s are FAA Compliant for in-flight use. They come with auto-adjusting voltage power supplies that let you plug into any outlet in the world with just an added plug adapter.

