

Clineb Basic

Aerosol delivery systems



CLINEB BASIC is a professional piston compressor nebulizer for aerosoltherapy. Indicated for intensive treatments in clinics or hospitals, its design includes a large compartment for convenient storage of accessories and medications.

AVAILABLE CONFIGURATIONS

- RE 300551 ————— Clineb Basic
- 27731 ————— Stand on castors (optional)

INCLUDED ACCESSORIES

HI-FLO Jet Nebulizer
 Adult Masks with elastic strings
 Child Masks with elastic strings
 Mouthpiece
 Nosepiece (non invasive)
 Air Tube
 Air Filter

FEATURES

Compressor ————— Oilless and maintenance-free piston pump
 Power Feeding ————— 230 V-50 Hz
 Power consumption ————— 184 VA
 Max Pressure ————— 2,5 bar 36 psi 250 Kpa
 Operating Pressure ————— 1,3 bar 18,9 psi 130 Kpa
 Max Air Flow ————— 16 L/min
 Operating Air Flow ————— 5,2 L/min
 Neb-rate with 4 ml solution ————— 0,5 ml/min
 Mass Median Aerodynamic Diameter — 3,25 μ m (MMAD)
 Respirable Fraction (< 5 μ m) ————— > 70 %
 Noise Level (EN 13544-1 Standard) — 57 dB
 Weight ————— 2,25 Kg
 Size (Main Unit) ————— cm 23x19x16
 Warranty ————— 2 years
 Place of Manufacturing ————— Italy