

Clineb Pro

Aerosol delivery systems



MAIN APPLICATIONS

Professional use

Intensive use

Room / ward

Professional piston compressor nebulizer for aerosoltherapy, indicated for intensive treatments in clinics or hospitals. Equipped with new breath enhanced jet nebulizer HI-4 with valve system reducing wasted volume and with 4 different speed settings. Position of the mask is adjustable by rotation allowing nebulization to bedridden patients or children.

AVAILABLE CONFIGURATIONS

● RE 300560 ————— CLINEB PRO

INCLUDED ACCESSORIES

1 x HI-4 Nebulizer
 1 x Mouthpiece
 1 x Nosepiece (non invasive)
 1 x Adult mask with elastic strings
 1 x Child mask with elastic strings
 1 x Air tube
 1 x Air Filter

FEATURES

Compressor ————— Oilless and maintenance-free piston pump
 Power feeding ————— 230V / 50Hz
 Max pressure ————— 3.0 bar 43.5 psi 300 kPa
 Operating pressure ————— 0.95 bar 13.7 psi 95 kPa
 Max air flow ————— 16 l/min
 Operating air flow ————— 8 l/min
 Neb-rate with 2ml NaCl 0.9% solution — position A (closed) 0.40 ml / pos. B 0.60 ml / pos. C 0.70 ml / pos. D 0.80 ml
 MMAD (µm) EN 13544-1 ————— position A (closed) 3.32 / pos. B 4.07 / pos. C 4.13 / pos. D 4.18
 Duty cycle ————— Non-stop operation
 Noise level ————— 60 dBA
 Weight ————— 2.2 kg
 Size (LxWxH) ————— 23 x 19 x 25 cm
 Warranty ————— 2 years