



GIMA

FAMILY ULTRASONIC NEBULIZER

Code: 28106
Category: Domestic use ultrasound nebulizers
Unit of sale: 1 pc.
Minimum order: 1
Type: Medical device
Class: II A
NSIS: 237511
CND: Z12040210
EAN13: 8023279281064



Description: ULTRASONIC NEBULIZERS

Thanks to the innovative electronic technologies these easy-to-use ultrasonic nebulizers provide a quiet, fast and portable operation with optimal particle size.

- **Portable:** small and lightweight, it allows patients to perform treatments anytime and anywhere.
- **Optimal particle size:** particle size determines the extent to which the medication can be inhaled by the patient into the deep lung tissue. Approx 90% of particles generated by Family or Pingoo have the optimal 1-5 μ (microns) size.
- **Fast treatment:** with its high-efficiency design, it allows a quick treatment time, most within 6 minutes, with limited power consumption.
- **Low noise level:** treatments are not only fast but are virtually silent due to its advanced ultrasonic technology.
- **Rechargeable battery:** optional
- **8 languages user manual:** GB, FR, IT, DE, ES, PT, GR, SA

Technical Specifications:

- **Operating voltage:** 220-240 V 50/60 Hz
- **Power:** 10 W
- **Ultrasonic frequency:** 2.5 MHz
- **Nebulization time:** 0.2-0.7 mL/min depending on medication density
- **Medication cups capacity:** max 8 mL
- **Automatic switch-off**
- **Particle size:** 5 μ (microns) average
- **Three programmable flow levels**
- **Noise level:** <40 dB

Standard accessories:

- **Adult facemask**
- **Pediatric facemask**
- **Nasal prong**
- **Mouthpiece**
- **Connection tube**
- **Medication cups**
- **Transformer (220-240V-12V DC)**
- **Carrying bag**
- **User manual**