

Miko

Aerosol delivery systems



MIKO is a piston compressor nebulizer for aerosoltherapy. Designed for home-care with compact size and carrying bag for easy transport of unit, accessories and medications. Recently upgraded for non-stop operation in case intensive treatments are required.

AVAILABLE CONFIGURATIONS

- RE 300600/03 ————— Miko
- RE 300600/12 ————— Miko Basic

INCLUDED ACCESSORIES

HI-FLO Jet Nebulizer
 Adult Masks with elastic strings
 Child Masks with elastic strings
 Mouthpiece
 Nosepiece (non invasive)
 Air Tube
 Air Filter
 Carrying Bag (Miko model only)

FEATURES

Compressor	Oilless and maintenance-free piston pump
Power Feeding	230 V-50 Hz
Power consumption	170 VA
Max Pressure	2,2 bar 32 psi 220 kPa
Operating Pressure	0,95 bar 12 psi 80 kPa
Max Air Flow	12 L/min
Operating Air Flow	4,8 L/min
Neb-rate with 4 ml solution	0,35 ml/min
Mass Median Aerodynamic Diameter	2,44 μ m (MMAD)
Respirable Fraction (< 5 μ m)	> 80 %
Noise Level (EN 13544-1 Standard)	56 dB
Weight	1,45
Size (Main Unit)	cm 14x17x11
Warranty	3 years
Place of Manufacturing	Italy