

# Mobile

## Aerosol delivery systems



MOBILE is a miniaturized micropiston compressor system for aerosol drug delivery featuring integrated "VRS" (Vibration Reduction System) and "RCC" (Rapid Charge Control) with electronic management of rapid battery charge.

It's a multipurpose device for home-care and hospital use, for both adults and children, and easy to carry anywhere.

The commercial rechargeable batteries can be purchased anywhere while the battery charger is integrated in the unit.

### AVAILABLE CONFIGURATIONS

- RE 300700 ————— Mobile

### INCLUDED ACCESSORIES

- HI-FLO Jet Nebulizer
- Adult Masks with elastic strings
- Child Masks with elastic strings
- Mouthpiece
- Nosepiece (non invasive)
- Air Tube
- Air Filter
- Carrying Bag

### FEATURES

Compressor	Oiless and maintenance-free piston pump
Power Feeding	100-240 V 50/60 Hz
	da 3,5 a 7 ore a seconda delle batteria usate
	40-90 min ( con batterie da 1600 a 2500 mAh
	8 x pile alcaline AA tipo 1,2 V
Power consumption	38,5 VA
Max Pressure	2 bar 29 psi 200 Kpa
Operating Pressure	0,8 bar 11,6 psi 80 Kpa
Max Air Flow	12 L/min
Operating Air Flow	4 L/min
Neb-rate with 4 ml solution	0,35 ml/min
Mass Median Aerodynamic Diameter	3 $\mu$ m (MMAD)
Duty Cycle	30 min ON- 30 min OFF
Respirable Fraction (< 5 $\mu$ m)	> 75 %
Noise Level (EN 13544-1 Standard)	50 dB
Weight	0,56 Kg
	12 V
Size (Main Unit)	cm 14x7x9
Warranty	5 years
Place of Manufacturing	Italy